

TPLO Saw Blade

Preoperative angle of the Tibial Plateau																										
	15°	16°	17°	18°	19°	20°	21°	22°	23°	24°	25°	26°	27°	28°	29°	30°	31°	32°	33°	34°	35°	36°	37°	38°	39°	40°
8 mm	1.4	1.5	1.7	1.8	1.9	2.1	2.2	2.4	2.5	2.6	2.8	2.9	3	3.2	3.3	3.5	3.6	3.7	3.9	4	4.1	4.3	4.4	4.5	4.7	4.8
12 mm	2.0	2.2	2.4	2.6	2.9	3.1	3.3	3.5	3.7	3.9	4.1	4.3	4.5	4.7	4.9	5.1	5.3	5.5	5.7	5.9	6.1	6.3	6.4	6.6	6.8	7
15 mm	2.6	2.8	3.1	3.3	3.6	3.8	4.1	4.3	4.6	4.9	5.1	5.4	5.6	5.9	6.1	6.4	6.6	6.9	7.1	7.4	7.6	7.9	8.1	8.4	8.6	8.8
18 mm	3.1	3.4	3.7	4	4.3	4.6	4.9	5.2	5.5	5.8	6.1	6.5	6.8	7.1	7.4	7.7	8	8.3	8.6	8.9	9.2	9.5	9.8	10.1	10.3	10.6
21 mm	3.6	4	4.3	4.7	5	5.4	5.8	6.1	6.5	6.8	7.2	7.5	7.9	8.3	8.6	9	9.3	9.7	10	10.4	10.7	11.1	11.4	11.8	12.1	12.4
24 mm	4.1	4.5	5	5.4	5.8	6.2	6.6	7	7.4	7.8	8.2	8.6	9	9.5	9.9	10.3	10.7	11.1	11.5	11.9	12.3	12.7	13.1	13.5	13.9	14.3
27 mm	4.7	5.1	5.6	6	6.5	7	7.4	7.9	8.4	8.8	9.3	9.7	10.2	10.6	11.1	11.6	12	12.5	12.9	13.4	13.8	14.3	14.7	15.2	15.6	16.1
30 mm	5.2	5.7	6.2	6.7	7.2	7.8	8.3	8.8	9.3	9.8	10.3	10.8	11.3	11.8	12.3	12.9	13.4	13.9	14.4	14.9	15.4	15.9	16.4	16.9	17.4	17.9

Rotation (mm) to provide the resulting 5° TPA